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Tecido Humano Doado

**READ BEFORE USING****DONATED HUMAN TISSUE****CAUTION: ALL ALLOGRAFTS ARE FOR SINGLE PATIENT USE ONLY****PASSED USP <71> STERILITY TESTS****THIS TISSUE WAS RECOVERED FROM A DECEASED DONOR FROM WHOM LEGAL AUTHORIZATION OR CONSENT HAS BEEN OBTAINED. THIS RECOVERY WAS PERFORMED USING ASEPTIC TECHNIQUES. PROCESSING AND PACKAGING WERE PERFORMED UNDER ASEPTIC CONDITIONS. TERMINAL STERILIZATION AGENTS WERE NOT USED IN THE PROCESS.****Description and Indication for Use**

MUSCULOSKELETAL TRANSPLANT FOUNDATION (MTF) tissues are supplied in a variety of standard sized units designed for surgical use by qualified health care professionals (e.g., physicians, dentists, and/or podiatrists). Processed human bone and soft tissue have been used in a variety of surgical applications and in combination with prosthetic devices. The amount and size of allograft necessary for a surgical procedure is based upon an individual surgeon's preference and the size and type of defect. The description of the tissue, serial number, expiration date, product code, size and/or amount, and additional information are printed on the allograft container label.

## **Contraindications**

- Fever
- Uncontrolled diabetes
- Pregnancy
- Hypercalcemia
- Renal insufficiency
- History of or active Pott's disease
- Sepsis, or infection in or around the surgical site
- Inability to cooperate with and/or comprehend post-operative instructions

## **Cautions and Warnings**

### **ALL ALLOGRAFTS ARE FOR SINGLE PATIENT**

**USE ONLY.** Do not use portions of an allograft from one container on multiple patients. Do not sterilize. Do not freeze. Trace amounts of Gentamicin, Primaxin, Amphotericin B antibiotics, Penicillin G, Streptomycin Sulfate, Dexamethasone Sodium Phosphate, vitamins, and Dulbecco's Modified Eagle Medium (DMEM) may be present. Tissue is exposed to processing solutions that may contain detergents and alcohol. Trace amounts of processing solutions may remain. Caution should be exercised if the patient is allergic to any of these substances.

Dispose of excess or unused tissue in accordance with recognized procedures for discarding regulated medical waste materials.

This allograft must not be used under any of the following conditions:

- If the container seal is damaged, not intact, or has physical damage.
- If the container has any physical damage.
- If the container label or identifying bar code is severely damaged, not readable or is missing.
- If the freeze-dried allograft container has been allowed to freeze or has otherwise been damaged.
- If the freeze-dried allograft has been rehydrated for more than 24 hours.

- If the expiration date shown on the container label has passed.

### **Precautions**

Extensive medical screening procedures have been used in the selection of all tissue donors for MTF (please see Donor Screening and Testing). Transmission of infectious diseases, may occur despite careful donor selection and laboratory testing, including serology and nucleic acid testing (NAT). Bacterial infection at the site of grafting may occur. Bacterial infection at the site of grafting may occur.

### **Adverse Effects**

Possible adverse effects of using human tissues include but are not limited to:

- Infection of soft tissue and/or bone (osteomyelitis)
- Immune response of non-infectious cause, including fever
- Deformity of the bone at the site
- Incomplete bone ingrowth, delayed union or non-union
- Fracture of the newly formed bone
- Disease transmission or undesirable immune response

**Adverse outcomes attributable to the tissue must be promptly reported to your local representative.**

### **Processing**

Processing and packaging are performed under controlled aseptic conditions in an ISO Class 4 environment.

- Tissue that is aseptically processed with no exposure to gamma radiation is labeled as follows: "Tissue is recovered and processed under aseptic conditions" and "Passes USP <71> Sterility Tests."

- Tissue that is aseptically processed and treated with low-dose gamma radiation is labeled as follows: "Tissue is recovered and processed under aseptic conditions. Treated with gamma radiation" and "Passes USP <71> Sterility Tests."

### **Donor Screening and Testing**

Prior to donation, the donor's medical/social history was screened for medical conditions or disease processes that would contraindicate the donation of tissues in accordance with current policies and procedures approved by the MTF's Medical Board of Trustees. Donor blood samples taken at the time of recovery were tested by a facility that is CLIA certified and registered with the FDA. The donor blood samples were tested for:

- Hepatitis B virus (HBV) surface antigen
- Hepatitis B virus core antibody
- Hepatitis C virus (HCV) antibody
- HIV-1/2 antibody
- HTLV-I/II antibody
- Syphilis
- HIV -1 NAT
- HCV NAT
- HBV NAT

All infectious disease tests results passed acceptability for screening. This allograft tissue has been determined to be suitable for transplantation.

The infectious disease test results, consent, current donor medical history interview, physical assessment, available relevant medical records to include previous medical history, laboratory test results, autopsy and coroner reports, if performed, and information obtained from any source or records which may pertain to donor suitability, have been evaluated by an MTF physician and are sufficient to indicate that donor suitability criteria current at the time of procurement, have been met. This tissue is suitable for transplantation.

The donor suitability criteria used to screen this donor are in compliance with the FDA regulations published in 21 CFR Part 1271 Human Cells, Tissues and Cellular and Tissue Based Products, as applicable. All procedures for donor screening, including laboratory testing, meet or exceed current standards established by the American Association of Tissue Banks.

### **Freeze-dried Bone**

All freeze-dried bone has been preserved using lyophilization (freeze-drying) to lower the total water content to 6% or less. All freeze-dried tissues are packaged in nested plastic trays.

### ***Preparation for use***

Preparation of the host bed is important for allograft incorporation. The host bed should be free of infection prior to grafting. Whenever possible, the allograft should be securely fixed to the host bone to aid in incorporation and to prevent displacement of the graft.

### ***Storage***

Store containers of freeze-dried tissue at ambient temperature. In order to maintain integrity of seal, do not freeze. It is the responsibility of the facility or clinician to maintain the tissue in the appropriate recommended storage conditions prior to use. If storage conditions or container seal have been compromised before intended use, the tissue should be discarded.

### ***Reconstitution/Rehydration procedure***

To obtain the best clinical results and prevent graft failure, the procedure and recommendations listed below should be followed.

The decision to rehydrate MTF freeze-dried bone prior to transplantation should be based upon the surgeon's preference. For tissues that are to

be cut, shaped, drilled or used for weight bearing purposes, excessive force should not be applied to the lyophilized bone during manipulation or upon implantation.

### ***Recommended instructions for handling freeze-dried tissue***

- Allograft tissue should be maintained in an aseptic environment at all times to prevent the possibility of contamination.
- It is common surgical practice to rehydrate freeze-dried tissue in an acceptable sterile irrigant (i.e. normal saline or Lactated Ringers Solution). Antibiotics may be used with the irrigant according to surgeon preference.
- Patient sensitivity to antibiotics used to rehydrate allograft tissues should be checked prior to use. Concentration of antibiotic solutions should be less than normally indicated for I.V. administration.
- Use new solutions for each allograft.
- Sufficient solution should be prepared to completely cover the tissue.
- Tissues should be implanted within 24 hours of opening the final tissue container as long as the allograft tissue is maintained in an aseptic environment.

### ***Opening of Package***

Note: The inner and outer tray components are sterilized. Use standard aseptic/sterile technique to open the package and make ready for use.

1. Peel back lid of outer tray. Note: Once the outer tray is opened, allograft should be used promptly. Inner tray, alone, although a sterile barrier, is not intended for storage of allograft, as it may not provide an adequate moisture barrier.
2. Grasp the pull-tab on the lid of the inner tray to remove it from the outer tray and pass it into the sterile field.

3. Peel back lid of inner tray. Transfer tissues to a sterile container for reconstitution.
4. Completely immerse the tissue in the reconstitution solution.
5. Rinse each tissue thoroughly with sterile irrigant prior to transplantation.

### **Patient Record**

Tissue recipient records must be maintained by the consignee and hospital for the purpose of tracing tissue post-transplantation. This will allow MTF to facilitate the investigation of actual or suspected transmission of communicable disease, and take appropriate and timely corrective action. A TissueTrace® Tracking Form and peel-off stickers have been included with each package of tissue. The serial number and the tissue description have been preprinted on the peel off labels. Please record the patient ID, the name and address of the hospital, allograft tissue information (using the peel off stickers), and specify the implant site or procedure performed using the bone allograft on the Tissue Trace Form. Once completed, the bottom page of the form should be returned to the local allograft representative or provider. Copies of this information should be retained by the hospital for future reference.

### **REFERENCES:**

Current MTF policies and procedures are in compliance with current FDA, AATB and other regulatory requirements.

### **Definitions of Label Symbols**

See IFU

Do Not Reuse

Facilitated through DePuy Synthes

Processed and distributed by: Musculoskeletal Transplant Foundation

125 May Street, Edison, NJ 08837, USA

1232 Mid-Valley Drive, Jessup, PA 18434, USA

+1.732.661.0202

All recovery, processing and distribution costs were paid for by MTF, a non-profit organization.

MTF Tissue forms and products are protected by one or more issued or licensed United States patents. A list of patents on available tissues and related technologies may be found on the MTF web site [www.mtf.org](http://www.mtf.org).

MTF Musculoskeletal Transplant Foundation® is a registered trademark of Musculoskeletal Transplant Foundation., Edison, NJ 08837 USA.

©2018 Musculoskeletal Transplant Foundation.

CTO 100024

FOR USE OUTSIDE OF THE UNITED STATES ONLY

*Please note: Human tissue for transplantation shall not be offered, distributed or dispensed for Veterinary Use*

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